

Changing Operation Settings on Thermocoagulator

1. Remove probe and insert battery
2. Hold button for 14 to 15 seconds until it flashes
3. Let go of button, wait another 12 to 13 seconds (it will stop flashing and then start flashing)
4. Wait until after the second flashing.
5. The unit will display a LED setting control option (see below).
6. Press the button to cycle through the options until you have the setting control option that you want.
7. Wait until it flashes and completely stops flashing.
8. The unit is ready for use. It will save the last saved settings for when it is turned on next.

Control Options for HTU

Time Temp	20sec	30sec	40sec	60sec
80°C	0 	1 	2 	3
90°C	4 	5 	6 	7
100°C	8 	9 	10 	11
120°C	12 	13 	14 	15

To Change the Option Setting, hold the ON switch down until the four Blue LED's begin to flash (~13sec), then release the switch to enter Option Selection Mode. (A probe is not required.) The first Mode Option is to change the sound setting. (If the switch is not pressed for ~6sec, then the four Blue LED's begin to flash again to enter the second Mode Option.)

Three sound options are available, normal sound during cycle (right LED On), no sound during cycle (no LED On), and an extra long sound at the end of the cycle (mid right LED On). Pressing the switch advances the sound options through the three choices. Continue to press and release the switch until the desired option is selected. No press of the switch for ~6sec, stores the last choice, and will advance the Selection to the second Mode Option, as indicated by Flashing of the four Blue LEDs.

The second Mode Option changes the Treatment time and temperature according to the Matrix above. Pressing and releasing the switch advances the option through the sixteen choices, (zero follows 16). Continue to press and release the switch until the desired option is selected, as indicated by the LEDs matching the Matrix. No press of the switch for ~6sec, stores this last choice, and will terminate the Option Selection Mode, as indicated by Flashing of the four Blue LEDs.